

PSG11-pregnancy specific beta-1-glycoprotein 11 Gene

Reference: PTXBC020711

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PSBG-11; PSBG-13; PSG13; PSG14; pregnancy-specific beta-1-glycoprotein 11; PS-beta-G-11; PS-beta-G-13; pregnancy-specific beta-1-glycoprotein 13; pregnancy specific beta-1-glycoprotein 11		

For research use only.

GENE DESCRIPTION	pregnancy specific beta-1-glycoprotein 11
GENE ID	5680
ORIGIN SPECIES	Human
ACCESSIONS	BC020711
VECTOR	pDONR223
NCBI SYMBOL	PSG11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacgcagcagagatcatggggccctctcagccctccctgcacagagcacatcaatggaaggggctcctgcacagtgagactcccaagccctccatccagcagcaactaaacccaggaggccatggagactgtgatcttaacctgtaacctgagactccggacgc
aagctacctgtggatgaatggatcagagcctccctatgactcataggaatgcagctgtctgaaacacagcaccctcttctatttgggtgcacaaagtatactgcaggaccctatgaatgtgaaatatggaactcaggagtgccagccgagtgaccagtcaccctgaatctcccat
ggfccagacctccccaatgttccctcagtcacctcttactatcaggagagaacctcgactgtctctgcagactcaccaccagcagatcttggacaattaatgggaagtttcagctatcaggacaaaagctcttatccctcaaatctcacaagcataatgggctctatgct
tgctctgctgtaactcagccactggcgaggaaagctccacatcttgacaatcagagtcattgctctccaggattaggaacttttgccttaataatccaacgtag