



Reference: PTXBC020722

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	methylglutaconyl-CoA hydratase, mitochondrial; 3-methylglutaconyl-CoA hydratase; AU RNA binding protein/enoyl-CoA hydratase; AU RNA-binding protein/enoyl-Coenzyme A hydratase; AU-binding protein/Enoyl-CoA hydratase; AU-specific RNA-binding enoyl-CoA hydra		

GENE DESCRIPTION	AU RNA binding protein/enoyl-Coenzyme A hydratase
GENE ID	549
ORIGIN SPECIES	Human
ACCESSIONS	BC020722
VECTOR	pDONR223
NCBI SYMBOL	AUH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1