

SH2D1A-SH2 domain protein 1A Gene

Reference: PTXBC020732

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SAP/SH2D1A; DSHP; EBVS; IMD5; LYP; MTCP1; SAP; XLP; XLPD; XLPD1; SH2 domain-containing protein 1A; Duncan disease SH2-protein; SLAM associated protein/SH2 domain protein 1A; SLAM-associated protein; T cell signal transduction molecule SAP; signaling lymph		

For research use only.

GENE DESCRIPTION	SH2 domain protein 1A
GENE ID	4068
ORIGIN SPECIES	Human
ACCESSIONS	BC020732
VECTOR	pDONR223
NCBI SYMBOL	SH2D1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgcagtggtgtgtatcatggcaaaatcagcagggaacccggcgagaagctcctgctgccactgggctggatggcagctatttgcaggacagcagagcgtgccaggcgtgtactgctatgtgtgtgtatcacggttacattatacataccagtggtccagacagaaa
cagggtcttggagtgctgagacagcacctgggtacataaaagatatctccgaaaataaaaaatctcatttcagcatttcagaagccagatcaaggcattgtaatacctctgcagtatccagttgagaagaagtcctcagctagaagtagacacaaggtactacaggataagagaagatcctga
tgtctgcctgaaagcccatga