

CYP2C9-cytochrome P450, family 2, subfamily C, polypeptide 9 Gene

Reference: PTXBC020754

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPC9; CYP2C; CYP2C10; CYP11C9; P450IIC9; cytochrome P450 2C9; cytochrome P-450 S-mephenytoin 4-hydroxylase; cytochrome P-450MP; cytochrome P450 PB-1; cytochrome P450, family 2, subfamily C, polypeptide 9; flavoprotein-linked monooxygenase; microsomal mono		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 2, subfamily C, polypeptide 9
GENE ID	1559
ORIGIN SPECIES	Human
ACCESSIONS	BC020754
VECTOR	pDONR223
NCBI SYMBOL	CYP2C9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattctgtggtcctgtgctctgtctctcatgtttgcttctcttctactctggagacagagctctggagaggaaaaactccctctgccccactctctccagtgattggaatctctacagataggtattaaaggacatcagcaaatcttaaccaatctctcaaaaggtctatggccctg
tgttcaactctgtattttggcctgaaccatagtggtgctgcatggaatgaagcagtggaaggaagccctgattgatcttggagaggagtttctggaagaggcatitccactggctgaaagagctaacagaggatttggattgtttcagcaatggaagaaatggaaggatccggcgt
ttctccctcatgacgctgcggaattttgggatgggaagaggagcattgaggaccgtgttcaagaggagcccgctgccttggaggagttgagaaaaaccaagggtgggtga