

FGA-fibrinogen alpha chain Gene

Reference: PTXBC020764

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Fib2; fibrinogen alpha chain; fibrinogen, A alpha polypeptide		

For research use only.

GENE DESCRIPTION	fibrinogen alpha chain
GENE ID	2243
ORIGIN SPECIES	Human
ACCESSIONS	BC020764
VECTOR	pDONR223
NCBI SYMBOL	FGA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgtttccatgaggatcgtctgcctggtcctaagtggtggcagcatggactgcagatagtggtgaagtgacttctagctgaaggaggagcgtcgtggcccaaggttggaagacatcaatcgcctgcaaaagattcagactggccctctgctctggtgaagactggaac
tacaatgccctctggctgcaggatgaaagggtgattgatgaagtcacaaagatttacaacagaataaataagctcaaaaatcactattgaatcagaagaacaataaggattctcattcgtgaccactaataataatgaaatttgagaggcattttcctcagccaataaccgtgat
aatcctacaaccgaggtcagaggatctgagaagcagaatggaagtcctgaagcgaagtcataaaaaattacagcaacaatttactagtagcaggttacacaagaggagactccacattgaaagcaagactataaaatggcagatgagccggaagtgaaagccgatcatg
aagaaacatagcaccagagaggccatgtaaatctgcctgtcagaggtatccactctccttggggaagccttccctgtccccctagactaagttaa
```