

GNG2-guanine nucleotide binding protein (G protein), gamma 2 Gene

Reference: PTXBC020774

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-2; g gamma-I; guanine nucleotide binding protein (G protein), gamma 2; guanine nucleotide binding protein gamma 2; guanine nucleotide-binding protein G(I)/G(O) gamma-2 subunit; G protein subu		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma 2
GENE ID	54331
ORIGIN SPECIES	Human
ACCESSIONS	BC020774
VECTOR	pDONR223
NCBI SYMBOL	GNG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagcaacaaccgccagcatagcacaagccagggaagctggtagagcagcttaagatggaagccaatcgcagagataaaggtgtccaaggcagctgcagatttgatggcctactgtgaagcacatgccagggaagacccctcctgacccctgttcggcttcagaaaaccgitttagggagaagaagttttctgtgccatccttaa