

CGA-glycoprotein hormones, alpha polypeptide Gene

Reference: PTXBC020782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CG-ALPHA; FSHA; GPHA1; GPHA; LHA; TSHA; glycoprotein hormones alpha chain; FSH-alpha; LSH-alpha; TSH-alpha; anterior pituitary glycoprotein hormones common subunit alpha; choriogonadotropin alpha chain; chorionic gonadotrophin subunit alpha; chorionic gon		

For research use only.

GENE DESCRIPTION	glycoprotein hormones, alpha polypeptide
GENE ID	1081
ORIGIN SPECIES	Human
ACCESSIONS	BC020782
VECTOR	pDONR223
NCBI SYMBOL	CGA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattactacagaaatatgcagctatcttctggtcacattgtcgggtttctgcatgttccattccgtcctgatgtgcaggattgccagaatgcacgtacaggaacccattcttctccagccgggtgcccaatacttcagtgcatgggctgctgtctctagagcatatccactccactaagggtccaagaagacgatgttgggtccaaaagaacgtcacctcagagtccttctgtgtagctaaatcatataacaggtcacagtaattgggggtttcaaaagtgagaaaccacacggcgtgccactgcagtactgtttattatcacaattcttaa