

NUMB-numb homolog (Drosophila) Gene

Reference: PTXBC020788

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NUMB, endocytic adaptor protein; numb homolog; h-Numb; protein numb homolog; C14orf41; S171; c14_5527		

For research use only.

GENE DESCRIPTION	numb homolog (Drosophila)
GENE ID	8650
ORIGIN SPECIES	Human
ACCESSIONS	BC020788
VECTOR	pDONR223
NCBI SYMBOL	NUMB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaattacggcaaaagttaggagaaagaggtgtttatgtccagagggcagtcgacacatcagtgccagacagatgaagaaggcgttcgaccggaaatgtagcttcccgttaagtaccttgccatgtagaagttgatgaatcaagagggaatgcacatctgtgaagatgctgtataaaagattgaaagctactggaagaaagcagtttaagcagttctgtgggtctcagcagatggactcagagttgtggatgaaaaactaaggacatagttgaccagacgatagagaagttcttctgtgcccagacaggaactttgatagagcctttctacatagccgtgatggcaccactcgtcgtggtatctgctcgttcatggctgtcaaggacacaggtgaaaggttgagccatgcagtaggctgtgcttttcagcctgttttagagcgcaagcagaagcgggagaaaggaatgtggagtgtactgtacttttgatgctagtcggaccacttttacaagagaaggatcattcgtgtcacacagccactgaacaagcagaagagagagatcatgaacaatgcaagatgccaagaaggtacaatctgtctacagtggtatttatggaatgcagttttcccaaggtatttttgaagctaccaggggtaaatttctcataatcaagagaaccttagcaataaa