



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SAA2-serum amyloid A2 Gene

Reference: PTXBC020795

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serum amyloid A-2 protein; serum amyloid A2		
For research use only.			

GENE DESCRIPTION	serum amyloid A2
GENE ID	6289
ORIGIN SPECIES	Human
ACCESSIONS	BC020795
VECTOR	pDONR223
NCBI SYMBOL	SAA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcttcctcagggccgtgtttctgctctctggttcctgagtcagcagccgaagctctttctgtcttggcagagcttttgatggggctgggacatgtggagagccctactctgcatgagagaagccaattacatcggtctcagacaaatactccatctctctgggggaactatgatgt
gccccaaaggagacctgggggtgctggtggccgcagaaagtcagcaatgccagagagaaatccagagactcacagcccatggttcggcaggagactcgctggccatcagctgcgcaataaaggggcaggagtggtgcagagaccccaactcacttcgacctgtgctgctcctgaga
aactatga