

PSMB1-proteasome (prosome, macropain) subunit, beta type, 1 Gene

Reference: PTXBC020807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HC5; PMSB1; PSC5; proteasome subunit beta type-1; macropain subunit C5; multicatalytic endopeptidase complex subunit C5; proteasome (prosome, macropain) subunit, beta type, 1; proteasome beta 1 subunit; proteasome component C5; proteasome gamma chain; pro		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, beta type, 1
GENE ID	5689
ORIGIN SPECIES	Human
ACCESSIONS	BC020807
VECTOR	pDONR223
NCBI SYMBOL	PSMB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtctctacagccatgtatcggctcctggcagagacttgggatggaaccgcacagagccggggcccttgcagctgcgattttccctacgttttcaacggaggtactactggcaattgctggagaagatttgcattgttctctgatactcattgagtgaagggtttcaat
tcatacggggatagcccaatgttacaattacagacaaaacagtcattggatgcagcgggttttcatggagactgtcttacgctgacaaagattattgaagcaagactaaagattgataagcattccaataaaggccatgactacggggcaattgctgcaatgctgtctacaatcctgt
attcaaggcgtctcttccatactatgtttacaacatcatcggtgacttgaagaaggaagggggctgtatcagctttgatccagtagggctctaccagagagactccttcaaggctggaggtcagcaagtgcctacagccctgctgacaaccaggttggttttaagaacatg
cagaatgtggagcatgttccgctgtccttggacagagccatgcggctgtgaaagatgtcttcttctgcgctgagagagatgtgtacactggggacgcactccggatctgcatagtaccaaagaggccatcaggagggaactgtttccttaaggaggactga