

TMEM68-transmembrane protein 68 Gene

Reference: PTXBC020835

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 68		
For research use only.			

GENE DESCRIPTION	transmembrane protein 68
GENE ID	137695
ORIGIN SPECIES	Human
ACCESSIONS	BC020835
VECTOR	pDONR223
NCBI SYMBOL	TMEM68
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatagacaaaatcaaacctgtggttaggacaggattctgtccctatatgattgtctgattcacatactcgaagaatgggttggtgaggagcagtgaggactattgaatttgcacactatctctgtgggttttacaccactaatacttttaacttcccttacttactatcttctctacc
ttactatttttttacacattataagagaagaagaatgtattgaagaagcctactctcataatttatgggatgtgcaaggaaaacagtggaactctgtgggatggacatgcagccgttggcatggttaagcaaggatactttcacctttgtgtgccattcatgtgtgctgcattggaactgttta
ccattccattttattgattga