

CLDN20-claudin 20 Gene

Reference: PTXBC020838

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	claudin-20; testis secretory sperm-binding protein Li 229n; claudin 20		
For research use only.			

GENE DESCRIPTION	claudin 20
GENE ID	49861
ORIGIN SPECIES	Human
ACCESSIONS	BC020838
VECTOR	pDONR223
NCBI SYMBOL	CLDN20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctcagcaggactccagctccttgccttcatctggccttatctggggtctctggagtcacagccactctgctcccaactggaagtgaaatgtggatgtggactccaacataacagccattgtacagctgcacgggctctggatggactgtacgtggtacagactgggatgt
tcagctgtgcctgaaacactccattctgtccctcccatccacgtgcaggctgcgagagccacatggtcctggcgtgtgttctgtctgtcttggggtctgcacttccacagttaggaatgaaatgtactcgttagggagggacagagaaaccaagagccatgctccttctgctgaggagt
ctgtttcatgtctgcaggaaatctctagttaatctgcagagtgtgtacacaaaggagatcatagcaaaacttctggatctgcaggttccagaaagcaacaacatgaacctggaggagctatctatcggattcattctgcaatgctgtgtttatctctggcatgtttctgcacctctgtataa
aaaggaatccagaagctagactcgacccacccacacagcagcctatctctaacacacagctcgagaacaattccacacaaatctgaaggattatgttaa