



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNASE6-ribonuclease, RNase A family, k6 Gene

Reference: PTXBC020848

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAD1; RNS6; RNasek6; ribonuclease K6; ribonuclease A D1; ribonuclease, RNase A family, k6; testicular tissue protein Li 166; ribonuclease A family member k6		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, k6
GENE ID	6039
ORIGIN SPECIES	Human
ACCESSIONS	BC020848
VECTOR	pDONR223
NCBI SYMBOL	RNASE6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgatgctttccctctttactctgctgctggtctatgggaccagtggtgtccactcatgtggtccaaagcgtctaccaaaggctcactggttgaaattcagcatatagccaaagcctctccaatgcaacagggcaatgagtggtcatcaacaattataccagcactgtaagcatcaaatactcttctgcatactcttcacaaatgctgctgctgctgtgatttctcagcatgtgctgcacaaatctgctgcacaaagcgtctcaaacagagcgtcaaaagcctgtcaacatcgtgactgactgcagactcaattcagggcatattcccaagtcgccgtatagtgctgctgcccaagtacaattcttcat