

SLC6A13-solute carrier family 6 (neurotransmitter transporter, GABA), member 13 Gene

Reference: PTXBC020867

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAT-2; GAT2; GAT3; sodium- and chloride-dependent GABA transporter 2; GABA transport protein; GABA transporter 2; solute carrier family 6 (neurotransmitter transporter), member 13; solute carrier family 6 (neurotransmitter transporter, GABA), member 13; s		

For research use only.

GENE DESCRIPTION	solute carrier family 6 (neurotransmitter transporter, GABA), member 13
GENE ID	6540
ORIGIN SPECIES	Human
ACCESSIONS	BC020867
VECTOR	pDONR223
NCBI SYMBOL	SLC6A13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatagcagggtctcaggcacaaccagtaatggagagacaaaaccagtgtatccagtcattggaaaagaaggagggaagatggcaccctggagcgggggcactggaacaacaagatggagtttgctgtcagtgctggggagatcattggcttaggcaacgtctggaggtttccct
atctctgctacaaaatgggggaggtgagatgagagccctgtgtccacccaccactcctggaaggagatacttcactctcctgcacttacggcccctctggggagtcctatag