

HIST4H4-histone cluster 4, H4 Gene

Reference: PTXBC020884

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4/p; histone H4; histone 4, H4; histone cluster 4 H4		

For research use only.

GENE DESCRIPTION	histone cluster 4, H4
GENE ID	121504
ORIGIN SPECIES	Human
ACCESSIONS	BC020884
VECTOR	pDONR223
NCBI SYMBOL	HIST4H4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcgaggtaagggtggcaagggtgggttaaggagggcgccaagcgccaccggaagggtgctggggacaatatccaaggcattacaagccggcgattcgccgacgtggggcggtcaagcgatttctgtctctacgaggagaccggggagtcct
caaagtctctggagaacgtgatccgtgacggtgactacacggagcagccaagcgcaagaccgtcacggccatggatgtgtgtacgcgctgaaacgccagggtcgccacctttatggttcggcggttga