

MGC24125-hypothetical protein MGC24125 Gene

Reference: PTXBC020886

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC728175		
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC24125
GENE ID	439935
ORIGIN SPECIES	Human
ACCESSIONS	BC020886
VECTOR	pDONR223
NCBI SYMBOL	MGC24125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaagaaagaagacgccacagacagcagcgagcagctaccacctccagagatgtgagtgacgctcttgacctctagccagctgcttcctgaatgcaatgcatgagcagggcaggcaaaactag