

TMEM140-transmembrane protein 140 Gene

Reference: PTXBC020942

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 140		
For research use only.			

GENE DESCRIPTION	transmembrane protein 140
GENE ID	55281
ORIGIN SPECIES	Human
ACCESSIONS	BC020942
VECTOR	pDONR223
NCBI SYMBOL	TMEM140
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggccaaggcctcagtgccgcgacgagctgtgtcatgagcatcatagtcctcgtgattgtgtgcatctcctgatgttatacgctcttctcgggaggtggcgaacctcaactgacctgccaacctgagaatcgctctataactctcctgtggaatgaggacaccagcaccc
tacagtgctaccagttccctgagctggaagccctgggggtgctcgggttgccctggccaggtggcgtgtacgggtcctcctcctcaccctcttggccccagcctcctcctagccccagtgcaacagtgatgagagagcgtggcggctggcagtgggcttcctgctgt
gtcctctgtcgtcgtggccggcggcctggtcctcctcctctctatgtgtggaagtgggtcaggctcctcctccggggcctgggtttctagctctgggcagcggccaggccttactcatcctcttgcctatagccatggcgtgtgtcctctgagggtgagaggcctgagagcaagcttgag
agctgctaa