



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM55B-transmembrane protein 55B Gene

Reference: PTXBC020947

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf9; type 1 phosphatidylinositol 4,5-bisphosphate 4-phosphatase; ptdIns-4,5-P(2) 4-phosphatase type I; ptdIns-4,5-P2 4-Ptase I; type 1 PtdIns-4,5-P2 4-Ptase; type I phosphatidylinositol-4,5-bisphosphate 4-phosphatase; transmembrane protein 55B		

For research use only.

GENE DESCRIPTION	transmembrane protein 55B
GENE ID	90809
ORIGIN SPECIES	Human
ACCESSIONS	BC020947
VECTOR	pDONR223
NCBI SYMBOL	TMEM55B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcagcatggagagcgttccccctgctgctgtagcccatcacgtgtggcgcggcgccgaacggtttatgtggggccggcgagggtgggctggggccgggggaggcctgacccctccgcacacacgtacggagccgattccccctgttccgaggggatccagcc
gtgttgcctggggaggaccacccccctattcaccttaactgacccggacatggggagtgccctatgatcacctgccgagtgctccaatctcatcaacgtgaaggcgaagatgcatcagcatgtagtcaaatgtgtgtctgcaatgaagcccccataagaatgacccccag
ggaaaaaatattgtcgtgcccctgtaactgtctcttattctgcaaatgtgacatcccaacggattgcatgccctggccctatgcaaaagaatcatcaacatggggcctgtgcatccggacacctgagtcgcgaaccccaacatgggtgtcagggttattctgacattgcaagaatac
ttttctgtagcagatgtcacagaccgacttttggctgtgtctcactgcgaagaaatgtcatctattggcgagatcacacgaataagatgtatctgctgtctgttggctgtcttggcctgactgactcaggctggccttggccttgcagagcatcagccgcgatgggc
actatgcacgcctggcaattctcatcctgtgtgctgtctgtctgttggcgccgctcttattggcgcttgaagtgctgaacgttccacgttccatgactctctcga

More info: Please contact us for any question: contact@proteogenix.fr