

SULF2-sulfatase 2 Gene

Reference: PTXBC020962

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSULF-2; extracellular sulfatase Sulf-2; sulfatase 2		

For research use only.

GENE DESCRIPTION	sulfatase 2
GENE ID	55959
ORIGIN SPECIES	Human
ACCESSIONS	BC020962
VECTOR	pDONR223
NCBI SYMBOL	SULF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaggcctcagtgcttcacccacgacaaccagcactggcagacggcgcttctggacactggggccttctgtgcctgcaccagcgccaataaacacgtactgggcatgaggaccatcaatgagactcacaattctcttctgtgaattgcaactggcttctagagtacttga
tctcaacacagacccctaccagctgatgaatgcagtgaaacacactggacaggggatgtctcaaccagctacagctacagctcatggagctgaggagctgcaagggttacaagcagtgtaaccccggaactcgaacatggacgtgggacttaagatggaggagctatgagcaataga
ggcagtttcagcgtcgaaatggccagaaatgaagagaccttcttcaaatcactgggacaactgtgggaagctgggaaggttaa