

YWHAG-tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, gamma polypeptide Gene

Reference: PTXBC020963

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	14-3-3GAMMA; PPP1R170; 14-3-3 protein gamma; 14-3-3 gamma; KCIP-1; protein kinase C inhibitor protein 1; protein phosphatase 1, regulatory subunit 170; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, gamma polypeptide; tyrosine 3-m		

For research use only.

GENE DESCRIPTION	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, gamma polypeptide
GENE ID	7532
ORIGIN SPECIES	Human
ACCESSIONS	BC020963
VECTOR	pDONR223
NCBI SYMBOL	YWHAG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggacgcgagcaactggtgcagaaagcccgctggccgagcaggcggagcgtacgacgacatggccggccatgaagaacgtgacagagctgaatgagccactgtcgaatgaggaacgaacaccttctgtgtggcctacaagaacgttggggggcagccgctcttc
ctggagcgtcatcagtagcattgagcagaagacatctgcagacggcaatgagaagaagattgagatggtccgtgctgacccgggagagatagagaaggagttggaggctgtgtgcccaggtgtgctgagcctgctggataactacatgatcaagaattgcagcagacccagtagcag
agcaaatgtgtctacctgaagatgaaaggggactactaccgctacgtgctgaaagtgccaccggagagaaaaggcgacggtgtgtggagctctctgagaaggcctacagcgaagccacgagatcagcaagagcagcatgcagccacccacccatccgattaggcctgctctt
aactactcgtcttctactatgagatccagaacgcccagagcaagcgtgcccactggccaagaccgcttcgacgacgccatgccgagcttgacacctcaacgaggaactcctacaaggactccacgctcatcagctcctccgcgacaacctcacgctctggacgagcagca
gcaggacgacgatggcggcgaaggcaacaattaa