



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF11-ring finger protein 11 Gene

Reference: PTXBC020964

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CGI-123; SID1669; RING finger protein 11

For research use only.

GENE DESCRIPTION	ring finger protein 11
GENE ID	26994
ORIGIN SPECIES	Human
ACCESSIONS	BC020964
VECTOR	pDONR223
NCBI SYMBOL	RNF11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaaactgctcaaatccccacctcggatgacatctccgtgcacagagctctagctccgaccgggctagcttggcgaggggacggagccggatcaggagcgccgccgcccatatcaggaaacaggtccagttcagctgtaccaccaacacttagccgactcggctagcaac
tcagctgactgaagagagcaaaattagatagatctcaagaagaattgtcttatacaacatctgcctaaaggaggtttatgacctgcgaagaagatgacatcagaagaataatcagaataaagatccggggagctgtgctgatctgtatgagtgactttgtatgggacaaacttgatttccgctgcatgcacatctatac