



Reference: PTXBC020965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	STT3A, catalytic subunit of the oligosaccharyltransferase complex; oligosaccharyl transferase subunit STT3A; STT3A, subunit of the oligosaccharyltransferase complex (catalytic); STT3A, cataylic subunit of the oligosaccharyltransferase complex; dolichyl-di		

For research use only.

GENE DESCRIPTION	STT3, subunit of the oligosaccharyltransferase complex, homolog A (S. cerevisiae)
GENE ID	3703
ORIGIN SPECIES	Human
ACCESSIONS	BC020965
VECTOR	pDONR223
NCBI SYMBOL	STT3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1