

## C12orf24-chromosome 12 open reading frame 24 Gene

Reference: PTXBC020967

|                |                                                                                                                                                       |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | uncharacterized protein C12orf24; C12orf24; HSU79274; protein FAM216A; protein predicted by clone 23733; family with sequence similarity 216 member A |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 12 open reading frame 24                                                               |
| <b>GENE ID</b>          | 29902                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC020967                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C12orf24                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtgtggacagctgtctccgcacacggctcgggtctcggcgcggagatgccggcagggtccgggtccgactggacggagcgtagctcttctgcagagccgcccgtgtggccgggacggagggtggcggcggcgatcagctggatactctgttaccagaattccaa  
agggttctgatagaatcaaaagatggatacaaaatgaactcacatagctaaagctgcaagagttatggaaaactccccaaatcaacaatccacctctctaaatcaatgatggaggcgtcttttcaagcatccagacctcaccagggccagaagcgttacctgtgcagcattgctaaaac  
tataatgcaactatctgaagatgttaagaggcagctacatgcagctactcagcagctcacaaaagccagggtgctctcactcatcacagaagccgcttagctcccgttactcacagaacagcattaccctgctactacatggcgacatcaactggagagaggactcgggtct  
tctgatatgcagctgcctgcactgaaatgctacacagcattccctttggcggccagtgagaacaaagaggataaaaactggatgcatctaaaacaagatgtaagtcactgaagattttagaagccaaggaaactgttcacaaacagtttctcagatgattctgaatcacac  
atgagtgaagaaaaaaggagaagatttactgaataattttatgcaatcaatgcaattgaagaacaggagaacatctgatgttaacttga