

PPP1R14A-protein phosphatase 1, regulatory (inhibitor) subunit 14A Gene

Reference: PTXBC021089

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPI-17; CPI17; PPP1INL; protein phosphatase 1 regulatory subunit 14A; 17 kDa PKC-potentiated inhibitory protein of PP1; 17-KDa protein; PKC-potentiated inhibitory protein of PP1; protein kinase C-potentiated inhibitor protein of 17 kDa; protein phosphatase		

For research use only.

GENE DESCRIPTION	protein phosphatase 1, regulatory (inhibitor) subunit 14A
GENE ID	94274
ORIGIN SPECIES	Human
ACCESSIONS	BC021089
VECTOR	pDONR223
NCBI SYMBOL	PPP1R14A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagctcagcggctgggcaagcgcgtgctgagcaagctgcagctccatgcgggcccgcgggcccagggggcagtcgggggggctgcagaagcggcagcgcgcgtcacctgaagtagaccggcgggagctgcagcggcggctggacgtggagaagtggatcgacg
ggcgctggaggagctgtaccggcgatggaggcagacatgccgatgagatcaacattgatgaattgttgagttagagagtgaagaggagagaagccggaaatccagggactcctgaagtcattgtggaaacctgtcgaggacttcacaggagctgtgcaaaagcttcaagg
cctccacaggcagccggcctccgccaagccctccacgacggcagcctcagccctccaggaccgggcccggactgtcacccctga