

TMEM208-transmembrane protein 208 Gene

Reference: PTXBC021109

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC171; transmembrane protein 208		

For research use only.

GENE DESCRIPTION	transmembrane protein 208
GENE ID	29100
ORIGIN SPECIES	Human
ACCESSIONS	BC021109
VECTOR	pDONR223
NCBI SYMBOL	TMEM208
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcccaagggcaagtgggcagagaggggaagaagcagatattgaagagaacagagagactctgaagtctacctgcggatcatactggggccaatgccattactgccttgtagcgttggtctctttactcatctgcctcatttgggcctggtggccctgggctttagtctggc
agtgtatggggccagctaccactctatgagctcgaatggcagcagcgttctctgaggatggggccctgatggatgggtggcatggacctcaacatggagcagggcagatgga