



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RBM4-RNA binding motif protein 4 Gene

Reference: PTXBC021120

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

LARK; RBM4A; ZCCHC21; ZCRB3A; RNA-binding protein 4; RNA-binding motif protein 4a; lark homolog; transcriptional coactivator CoAZ; zinc finger CCHC-type and RNA binding motif 3A; RNA binding motif protein 4

For research use only.

GENE DESCRIPTION	RNA binding motif protein 4
GENE ID	5936
ORIGIN SPECIES	Human
ACCESSIONS	BC021120
VECTOR	pDONR223
NCBI SYMBOL	RBM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtagaagctgttcaatcggaacactgccccggggagctacagagcagggaattcgcctcactctcgagcagtagtgggaaggtgcgtgaatgtgacatcaltgaagaattacggcttftgcacatagaagacaagacggcagctgaggatgccatcacgaacctgcacattacaagcttc
atgggttgtaacatcaacgtggagccgacgaagataacagcaaaactcaaaaggttcgatctggcgcaacatcagctccactcgacaaataaggactcgtgagccaaattggagggaatgtgcgtctgcacgaatgtgacatgctgaaagattatgcctctgacatcaatgaagcggcgg
cagagatgcagttggggcccatcagctctgtgataacacagcagttcaagttgaaccacctcttggtagagaggtctgtgacaaacagctgtgtgcagaaaatgtgtgtgaataaataacacagctgtgcctgctgctgtatgaagaataata