



Reference: PTXBC021121

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine-protein kinase WNK1 2; serine/threonine-protein kinase WNK1 1; serine/threonine-protein kinase WNK1; HSN2; HSN2; KDP; PPP1R167; PRKWKN1; PSK; p65; WNK lysine deficient protein kinase 1 isoform; erythrocyte 65 kDa protein; prostate-derive		

GENE DESCRIPTION	WNK lysine deficient protein kinase 1
GENE ID	65125
ORIGIN SPECIES	Human
ACCESSIONS	BC021121
VECTOR	pDONR223
NCBI SYMBOL	WNK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]