

## SMAP2-small ArfGAP2 Gene

Reference: PTXBC021133

|                |                                                                                                                       |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SMAP1L; stromal membrane-associated protein 2; stromal membrane-associated GTPase-activating protein 2; small ArfGAP2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | small ArfGAP2                                                                                    |
| <b>GENE ID</b>          | 64744                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC021133                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SMAP2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgacaggcaagtcggtgaaggacgtggatcggtaaccaggctgtcctggccaacctgtgctggaggagataacaagtttgtgcagattgccagtcctaaaggccgcgatggccctcttggaaattggtgtgtcatctgcattgatgtgctggaatccacaggaatctgggggtgca  
catatccagggttaagtcagtttaacctcgaccagtggactcaagaacagattcagtcgcatgagagatgggaaatggaaggcacaaccgactttatgaagcctatctctgagaccttccggcgacctcagatagaccagctgttgaaggatttattcgagacaaatagagaaga  
aatacatggaccgaagtcctggacatcaatgctttaggaaagaaaaagatgacaagtggaagaggagcgaaccagttccagaaaaaaatggaaacctgtgttttggagaagtgaaaatgccacagaaaaaagagaccacagctacctcggaaaagctccccgaaatccaca  
gcgcctgtcatggattgttggcccttgatgctcctgtggtcctgtcctcattgcaaatagtaagaccagcaataacctagagaaggatttagatctgttggcctctgttccatcccccttcttccggttccagaaaggtttaggttccatgccaaactgcaggagtgccggtctgttctgaaaa  
tctgaacctgtttccggagccaggagcaaatcagaagaataggaagaaacagctctctaaagactcattcttcatctgtatggatcccagagcctcaaatgctactcaagcaatgttcatggctcccgtcagatggcatatcccacagctaccccagctccccgggttacacct  
cctaacagcataatggggagcatgatgctccaccagtaggcaggttctcagccaggagcttctggatgttgcctccatggccatgctgcagcctatatgggtggcatgcaggcatcaatgatgggtgtgccgaatggaatgatgaccaccagcaggtggtctacatggcagg  
catggcagctatgccccagactgtgtatgggtccagccagctcagcagctgcaatggaaccttactcagatgacccagcagatggctgggatgaacttctatggagccaatggcatgatgaactatggacagtcgaatgagtgccggaatggacaggcagcaaatcagactctcagtc  
tcagatgtggaataa