

## DDI1-DDI1, DNA-damage inducible 1, homolog 1 (*S. cerevisiae*) Gene

Reference: PTXBC021172

|                |                                                                                                                                   |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | DDI1, DNA-damage inducible 1, homolog 1; protein DDI1 homolog 1; DNA-damage inducible protein 1; DNA damage inducible 1 homolog 1 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | DDI1, DNA-damage inducible 1, homolog 1 ( <i>S. cerevisiae</i> )                                  |
| <b>GENE ID</b>          | 414301                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC021172                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | DDI1                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtgtacacccgtgtactgctgctggaggacgtctccgaggtcacctctctccaggtcagccccgattttgagctccgaaactcaaggtcctctgctgaagcggagtcagagtcctccgctgaagagatccagatccacatggagcgaactctcatcaggaccactgttcct  
gggctctacggcctcaaatggtcgatctgtgttttactgcagaaggacaatgtgggacctcgggctccagggcgtgccccgaaccagcctcgtgtagacttcagtggtcattgctggtgctgggagctccagctccgctccacagcacctggacagcagcagcgcacaccc  
gctgcccagcgggtcacagggtgtgctgctcaggagagaaggtggccgctgcaaggtctgggctgccccgctgacccagcagatgctgctcctcaacccccacgactgtccctgtcctcaaggaacgcaacccctcctggcggaagcctgtcagcggaaagccttgagaccttt  
ctcaggtgctgtgagcagcacaagggaagggccttgagagagcaagagagccttctctctacacagccgacccactggatcgggaagctcagggccaaatagaagaggaatcggcagcaaaacattgaagaaacatgaatatagcagatagaagggccccgagagtttt  
gaacaagtgcagctctcatttaactgcaaatggaatgggcatcctttgaaggctttgtgactcgggccccagatgaccattatgagccaggtgtgctgagcagatgtaacatcatgaggtgtggtgacccagcgtgggtgggtgtgctaaaggagtgggcacacagagaattat  
tgcccggttcacatcagctcagattcaaatgaaggtgatttctacagtgctcttctccatacttgaggatcaaccatggatgctcttagcctagatgctccggagacatcaatgttccatcgattgaagaaaaatgtgctggtcatcggcaccactggcagcagacttatcttctctg  
agggagaggtgccccttatgctctaggatggttaagtgggcaagatgagttctcgacaaggaaattacacattcagtcagtgattcaggacgaaaagaacattaa