

ZCCHC13-zinc finger, CCHC domain containing 13 Gene

Reference: PTXBC021176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNBP2; ZNF9L; zinc finger CCHC domain-containing protein 13; zinc finger, CCHC domain containing 13; zinc finger CCHC-type containing 13		

For research use only.

GENE DESCRIPTION	zinc finger, CCHC domain containing 13
GENE ID	389874
ORIGIN SPECIES	Human
ACCESSIONS	BC021176
VECTOR	pDONR223
NCBI SYMBOL	ZCCHC13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcagtaaggattcttcggtgtgacactctggccattgggctcgggatgccctagaggagagctggaggcggaagaggtggaggccatggcagaggttctcaatgtggtccaccacccctatcttacacctgttactgtgtggtgagtcgggtcgtaatgctaagaactgtgt
ccttctcggaaacatctgtacaaactgtgggagaaagcggccacatcgccaaagactgtaaggatcctaaacgagagagacgccaacactgttatacctgcccagactaggacatctggctcgtgactgtgatcgtcagaaagagcagaaatgctactcttgcggcaacttgggcacatt
cagaaagactgcgccagggtcaaggtttaccgatgcggcgagattggccacgtggccatcaattgcagcaaggcgaggccaggtcaactgctaccgtcggcgaatccgacatctagccaaggaaatgtcccagtgta