

TCTEX1D2-Tctex1 domain containing 2 Gene

Reference: PTXBC021177

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tctex1 domain-containing protein 2; Tctex1 domain containing 2		
For research use only.			

GENE DESCRIPTION	Tctex1 domain containing 2
GENE ID	255758
ORIGIN SPECIES	Human
ACCESSIONS	BC021177
VECTOR	pDONR223
NCBI SYMBOL	TCTEX1D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacgtccatcggagtgctcttcggtggcgacgggtgcctgagcctgagaagaacgcaggggagcccgagaacacctattctcggcctgttttcagcagaggttcaggccctctgtggttaagactgtatccatgctgtgctcaaggaggaggactggcaaatgctgaat
attctccagaagaatgctcagcttacaacattatcagaaaacattaaagataaattaaaagaatgggatttgaccgatacaaaatgggtgcaagtagtgattggagaacaagaggtgaaggagatttcattgcttctcgtgttctgggatgctgacactgacaactatactcatg
atgtttcatgaatgacagttttatctgcgttgtagcagcatttgctgtttctactactga