

C16orf78-chromosome 16 open reading frame 78 Gene

Reference: PTXBC021181

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C16orf78; chromosome 16 open reading frame 78		
For research use only.			

GENE DESCRIPTION	chromosome 16 open reading frame 78
GENE ID	123970
ORIGIN SPECIES	Human
ACCESSIONS	BC021181
VECTOR	pDONR223
NCBI SYMBOL	C16orf78
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagagcaacaatggacctgaaggatttaatgccacaagaggaaatcatgtggaagactgctgaagataggcgcatgtctgacctcactgtgtgctggagtgctggagcggagggcaggggaagaagaacaagctcccgagaagcaaaagccaaagtggtagacagtc
cttaaacgaaataagaagaaggagaagaaggaagcaaaaggccctcatgacagcacggggagggaacccgagggacacggagacttccagcaggcccttaggaaagagattcagggaaggagccgcctctaccgaagcctctatggagtgaggcaaaaggggaacacacacag
catggtccctggcagctacatcaaggatggcccaagaatctgatacagacatcaaggatgcagtcgacccagagtgccactcagcggccaacccatccgacgacaaagcattgtcttagatcccatgttacaggagggtacaccttaacagccagaggccaaccttcataagagactgg
tccaacaagatgctgacatggttacgaacgcaagctaaagagcctcatggagaaaagcaccgaacctaaagatggaacatgaggatgttgaagccagaggaggtgctgagctgcccatacctgaggttatccaaggagaacattggaccttgctcaagttgtgcaaggatgcagg
aatgaatgtggatatccaccccatggtcgaagaggacatagatgctaaaaaggtgttcaccggaataccagcatggcctctaa