

DKFZP434K028-hypothetical LOC26070 Gene

Reference: PTXBC021187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC26070		

For research use only.

GENE DESCRIPTION	hypothetical LOC26070
GENE ID	26070
ORIGIN SPECIES	Human
ACCESSIONS	BC021187
VECTOR	pDONR223
NCBI SYMBOL	DKFZP434K028
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttcaggccttgacacagaaccacagaacatgctcaacaagcccccttaagcttaggggcactgacctacccctggggactagtaattttcaagatgagaaagcggagggcacagcgggtgccctgatacggccaagagaagcactgtgggttaagcaggagcctaggaacagacc
cagctgccctcgtcggcgaatgcacagaaccattacatctaccccgagaccaaccaagaggtacccacacccgtggcagccggggctgaaggcctggagataggagaggagggcccccatactccagcaccgaaggcaatgagagaactgaggccaggaaggtgcta
ggtcctgaagaggaaggtgggaaggcctgggcagtga