



Reference: PTXBC021193

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTC-479C5.17; KCC1; hKCC1; solute carrier family 12 member 4; electroneutral potassium-chloride cotransporter 1; erythroid K-Cl cotransporter 1; potassium/chloride cotransporter 1; solute carrier family 12 (potassium/chloride transporter), member 4; solut		

For research use only.

GENE DESCRIPTION	solute carrier family 12 (potassium/chloride transporters), member 4
GENE ID	6560
ORIGIN SPECIES	Human
ACCESSIONS	BC021193
VECTOR	pDONR223
NCBI SYMBOL	SLC12A4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]