



Reference: PTXBC021201

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JP-45; JP45; junctional sarcoplasmic reticulum protein 1; 2310032K21Rik; homolog of mouse skeletal muscle sarcoplasmic reticulum protein JP-45; junctional-face membrane protein of 45 kDa homolog; skeletal muscle sarcoplasmic reticulum protein JP-45		

For research use only.

GENE DESCRIPTION	junctional sarcoplasmic reticulum protein 1
GENE ID	126306
ORIGIN SPECIES	Human
ACCESSIONS	BC021201
VECTOR	pDONR223
NCBI SYMBOL	JSRP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]