



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

EIF4H-eukaryotic translation initiation factor 4H Gene

Reference: PTXBC021214

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WBSCR1; WSCR1; eIF-4H; eukaryotic translation initiation factor 4H; Williams-Beuren syndrome chromosome region 1		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 4H
GENE ID	7458
ORIGIN SPECIES	Human
ACCESSIONS	BC021214
VECTOR	pDONR223
NCBI SYMBOL	EIF4H
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1