

C7orf33-chromosome 7 open reading frame 33 Gene

Reference: PTXBC021251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf33; chromosome 7 open reading frame 33		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 33
GENE ID	202865
ORIGIN SPECIES	Human
ACCESSIONS	BC021251
VECTOR	pDONR223
NCBI SYMBOL	C7orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaagtggaggtcaagcctcagccttgaagagtgtccctggagactccaggcccccaatgtgaatgtgaagccctcctgccagtgaggcgaaggcgcggattgaccttcgctgagtgaggaggcagtcgctgtttgggtccacgttagggcggtccagggtcaatttaactgt
catatgccacgggaagacataagaagccaaacccacacaaaacatgaaccggggatggaatttattgtcctgtatcagctcccacaaatctgggtcccgctggcattttcttcaaggtcccacggatgccagagagcagtcagaatcaggccaggcaccaggatggcgtgtc
ctcagatccagttgtcggcacctgtctccagttacctagatctgctaactctcttacaagcctgggaggacagtgacaagttcatacctaaacgtacgtggcgtgaagtaagaaagctccagaattctgttgaggccacaaggattccagaactgacagcagttaa