

HMG20B-high-mobility group 20B Gene

Reference: PTXBC021585

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRAF25; BRAF35; HMGX2; HMGXB2; PP7706; SMARCE1r; SOXL; pp8857; SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1-related; BRCA2-associated factor 35; HMG box domain containing 2; HMG box-containing protein 20B;		

For research use only.

GENE DESCRIPTION	high-mobility group 20B
GENE ID	10362
ORIGIN SPECIES	Human
ACCESSIONS	BC021585
VECTOR	pDONR223
NCBI SYMBOL	HMG20B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggcgccgagtgaggcaagctgcagccaacggaaaagcagcggtacctggatgagggcagagagagaagcagcagctacatgaaggagctgcccgtaccagcagctgaagcctataagatgtgcacggagaagatccaggagaagaagatcaaaaaaagactc
gagctctgggctcatgaacactctcctgaatggacacaagggtggggaactgcgatggcttctccaccttcgatgttccatcttcaactgaagagttcttggacccaaaacaaagcgtgagggcggagcttcggcgcttgcggaagatgaatgtggccttcgaggagcagaacgcggtaactg
cagaggcacacgcagagcatgagcagcgcgcgagcgtctggagcaggagctggcgctggaggagcggaggacgctggcgctgcagcagcagctccaggccgtgcccagggcgtcacgcgagcttcgctcactgcgggtgcccgggcacgggcaaacgccacgctg
ggcactctggacttctacatggccggcttcacggagccatcgagcgcgacccgccagcagcagaagctcatctccgcatcaaggaaatctggcccaggtgcagcagcagcctgtga