

SLC16A4-solute carrier family 16, member 4 (monocarboxylic acid transporter 5) Gene

Reference: PTXBC021664

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MCT4; MCT5; monocarboxylate transporter 5; MCT 4; MCT 5; monocarboxylate transporter 4; solute carrier family 16 (monocarboxylic acid transporters), member 4; solute carrier family 16, member 4 (monocarboxylic acid transporter 5); solute carrier family 16		

For research use only.

GENE DESCRIPTION	solute carrier family 16, member 4 (monocarboxylic acid transporter 5)
GENE ID	9122
ORIGIN SPECIES	Human
ACCESSIONS	BC021664
VECTOR	pDONR223
NCBI SYMBOL	SLC16A4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgaagaggagggaaggtccaaccttacactaaacccctggatgggggatggggatggatgattgattcatttttctcctggatgaatgtgttgatggggatgaccaagacttttgaatttcttgggtctttcaagaagagtttgaaggcacctcagagcaaatgggtggattg
gatccatcatgcatctctctgttttggcagggtccctgggtgctattatttggacatactggagagaaaactaccatcttggggcttctgttactggtggatctgacagcagctgggccacaagtattcttttgggtgactatgggacttctaccgggttgggttctgttctt
ataccaagtggctgctggtgaactaccaatactcaaaaaacgattggctctttctacagctattgcccgttctgggatgggactgacttttcttggcaccctttacaaaattcctgatagatctgtatgactggacaggtatccctggagacgggtcagtcagattttctggatgggtgctgac
aaaactggattagaagatattaccacaagctttacctcctctcgcggcatcaactaactgcttgcctcttagccaccacattccactactatgacctacacatctgcttgcacatcttctggtggttacctggcattgatactgcctgactgggtgatctgttaggaattctacagtaa
acagggttttgggacttccagtttcttgggtggtgctgctcttctggaccacatagcaggtggttatatgattataccagacatacaatggctcttctactctctgcatatgctatctctctctcagtttcttctttttgtaccattggccgaaagatggaaaaacagctgacctga