

TSPAN2-tetraspanin 2 Gene

Reference: PTXBC021675

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET3; TSN2; TSPAN-2; tetraspanin-2; new EST tetraspan 3; tetraspan 2; tetraspan NET-3; tetraspan TM4SF; tetraspanin 2 isoform; tetraspanin 2		

For research use only.

GENE DESCRIPTION	tetraspanin 2
GENE ID	10100
ORIGIN SPECIES	Human
ACCESSIONS	BC021675
VECTOR	pDONR223
NCBI SYMBOL	TSPAN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcgttccgccccggcctgcggtgcatcaagtacgtgctgttggctcaacctgctcttctgctgctggatcggcgtcattgctttggactatggttcggttcggaggtgccataaaggagttatcatcagaggacaagtcgccagagtatttctatgtgggctgtatgtctggttgagccgggcccctgatgatgcccgtgggttcttctggatgctgcggagccatgcgggagtcgcaatgtgtgcttggatcatttttacctgctcctggtgataattgctgctgaaagtaaccactggagttattgcttttataggcaaggggtagctatccgacatgttcagaccatgtatgaagaggttacaatgattaccttaagacaggggaaaggcaatgggacactcatcacttccactcaacatttcagtgctgtgaaaaaagctccgaacaggtccaacctacatgccaaaggagcttctaggacacaagaattgcatcgaattgagaccataatcagtgtaagctccagctcattggaattgctgattggaattgcaggtctgacgatcttggcatgatattcagcatggtcctctgctgtgcgatacgaactcacgagatgtgatga