

39706-15 kDa selenoprotein Gene

Reference: PTXBC021697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	15 kDa selenoprotein; SEP15; selenoprotein F		

For research use only.

GENE DESCRIPTION	15 kDa selenoprotein
GENE ID	9403
ORIGIN SPECIES	Human
ACCESSIONS	BC021697
VECTOR	pDONR223
NCBI SYMBOL	39706
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggagctattcttgaaagttgtggatgaaattgggaaggtccctcaagtcgaagtgtgctgtggtcagaccctgtattaaagctttggacgacaaatgggaacattgctgaagaactgagcattctcaaatggaacacagacagtgtagaagaattcctgagtgaaaagttggaacg
catataa