

C14orf48-chromosome 14 open reading frame 48 Gene

Reference: PTXBC021728

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf48; c14_5713; long intergenic non-protein coding RNA 521		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 48
GENE ID	256369
ORIGIN SPECIES	Human
ACCESSIONS	BC021728
VECTOR	pDONR223
NCBI SYMBOL	C14orf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacacagccagagagctgaccaagcaatctggtgacaagggaagggaccttcaagggctgtggcagggaactctaccagctcaggctaagcagaagaagctcaagagagaagtcgagaagcacaagcttttgaagactatctgattaaggctctgagaaaatccccgagg
gctgcacgggatgggagccggagggaggtgctggggagggccacgggtgaagcactacgggaagctctcacagccagccagggacacgagaagcgctcaggccctctgccagatgatccaggctgtccacggagcctggagctctggagaggaccacagggctctcat
agcctcaagatccggctgtgtcagctgcagaagaagtgctaccgcaagcaggagcagtggtggcagctga