

EEF1DP3-eukaryotic translation elongation factor 1 delta pseudogene 3 Gene

Reference: PTXBC021729

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EF1G; GIG35; elongation factor 1-gamma; EF-1-gamma; PRO1608; eEF-1B gamma; pancreatic tumor-related protein; translation elongation factor eEF-1 gamma chain; eukaryotic translation elongation factor 1 gamma		

For research use only.

GENE DESCRIPTION	eukaryotic translation elongation factor 1 delta pseudogene 3
GENE ID	196549
ORIGIN SPECIES	Human
ACCESSIONS	BC021729
VECTOR	pDONR223
NCBI SYMBOL	EEF1DP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctacaaacttctagtaggtgagaagatctggtccacaagttcaatatggcgatgcagaaaggagattctacgaacagatgaacgggctgtggccggcgctccctccaggaggccagcatgatcctccatgatattgccagagccagagagaaatccgaaatccctggcc
ggagccttagggccaggggctctagcggccccagcggagaccacagcagctcgtctccggatcgcagctctggaagtggacaaccagagagacctggctgaacgtgctggagaagagctggccggccactgggccacagccccgagaccagcacatgtctccatgcg
ccaagtggagccccggccaagaagctagccacaccagcagagatga