

## ZNF248-zinc finger protein 248 Gene

Reference: PTXBC021819

|                |                                                                                                                              |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | bA162G10.3; zinc finger protein 248; KRAB protein domain; Kruppel-type zinc finger protein; bA162G10.3 (zinc finger protein) |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger protein 248                                                                           |
| <b>GENE ID</b>          | 57209                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC021819                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ZNF248                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaacaaatcccaggacaagtgtcattcaaggatgtatgtgtggacttcactcaggaagtggtatctgtggaccctgtcagaagattctacagagatgtgatcctggaaaattatagcaatctgtctcagtagggattgcattactaaaccagaagtgatctttaagatcgagcaa  
ggagaagagccctggatattagaaaaaggattcccaagccagtgccaccagaaaggaaatggaaagttgatgacgtgttagagagcagccaggaaaatgaagatgaccattttgggagcttctattccacaacaacaaacacattga