

RPL27-ribosomal protein L27 Gene

Reference: PTXBC021886

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L27; 60S ribosomal protein L27; ribosomal protein L27		
For research use only.			

GENE DESCRIPTION	ribosomal protein L27
GENE ID	6155
ORIGIN SPECIES	Human
ACCESSIONS	BC021886
VECTOR	pDONR223
NCBI SYMBOL	RPL27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcaaaacaaagcatctaaaaccgagtttctggaagaaccactgtctggctggagcgtactccggacgcaaaagctgtcatcgtgaagaacattgatgatggcacctcagatgccccacagccatgctctggtggtggaattgaccgctacccccgaaagtgcagctgcat
gggcaagaagaagatcgccaagagatcaaaataaatctttgtgaaagtgtataactacaatcacctaataatgccacaagggtactctgtggatccccctggacaaaactgtctcaataaggatgtcttcagagatcctcttaaacgcaaggccccgacgggagccaaggtaagt
tgaagagagatacaagacaggcaagaacaagtgtgtcttcagaaactgcggttttag