



Reference: PTXBC021965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B3GNT1; B3GN-T1; B3GNT6; BETA3GNTI; MDDGA13; iGAT; iGNT; beta-1,4-glucuronyltransferase 1; N-acetyllactosaminide beta-1,3-N-acetylglucosaminyltransferase; UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 1; UDP-GlcNAc:betaGal beta-1,3-N-acetylgl		

For research use only.

GENE DESCRIPTION	UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 1
GENE ID	11041
ORIGIN SPECIES	Human
ACCESSIONS	BC021965
VECTOR	pDONR223
NCBI SYMBOL	B3GNT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]