

C6orf64-chromosome 6 open reading frame 64 Gene

Reference: PTXBC022007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf64; SAYSVFN domain-containing protein 1; SAYSVFN motif domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 64
GENE ID	55776
ORIGIN SPECIES	Human
ACCESSIONS	BC022007
VECTOR	pDONR223
NCBI SYMBOL	C6orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacagcggttagctgagtttggcgggcggaacggcggttctggcgcccaaccctgctccagtcaggcgccacaacccaggagagaagcggaagcagcagcagctctaaaggcagccccaggctggctaaagcggttctggtaggaaacctaggccc
gcgagtgccccggccagcccgctagttcaggaagcggctcagccccaggcgacacatcagagacacatggaacacagccattctctgctgctgggaccagctcttctgaccaataacacctctgaaggtctctctggttgctcctgctgggactgttggaactg
gaattggcctggcctatgttctgctctgtctatgtggtgacgtcgggacacagggccctgaagagaagaagaggagagaagagcgctactctgtgttcaatccaggctgtgaagccatccaggccacccctgactgcagagcagttggagcgcgagttacagttgagaccct
ggcaggagatag