

LOC541473-FK506 binding protein 6, 36kDa pseudogene Gene

Reference: PTXBC022013

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FK506 binding protein 6, 36kDa pseudogene
GENE ID	541473
ORIGIN SPECIES	Human
ACCESSIONS	BC022013
VECTOR	pDONR223
NCBI SYMBOL	LOC541473
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggggaagcgcgttaaacaggagtcctggaaggggacgacgccccggccagtcctgtacgagcgggttaagtcagaggatgctggacatctcgggggaccggggcgtctgaaggacgtcatccgagaaggagctggagacctagtgccgcctgatgctcggtgctag
tgaaatactatggatacctggaacacttggacagacccttcgattctaattactttaggaaaactctcggcctaataatgaaacttggagaggatattacattgtgggcatggagctgggccttctgagcatgcagagaggagagctggccagatgcttcgtcttgggtaaactctcgactccaa
ggcccagcctccatctttacctcagagcctctga