

EXOSC1-exosome component 1 Gene

Reference: PTXBC022067

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-108; CSL4; Csl4p; SKI4; Ski4p; p13; exosome complex component CSL4; 3'-5' exoribonuclease CSL4 homolog; exosomal core protein CSL4; homolog of yeast exosomal core protein CSL4; exosome component 1		

For research use only.

GENE DESCRIPTION	exosome component 1
GENE ID	51013
ORIGIN SPECIES	Human
ACCESSIONS	BC022067
VECTOR	pDONR223
NCBI SYMBOL	EXOSC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccacctgtgagatactgcatccccggcgaaactgtgtgtaacttgaggaggggcagccggggcagcgccacacacccgggtacatcttttcgtcgttgcgggtgtctgatgaagagcagcgagaatggcggttcctcagtggtgtctgtagtggagaaacaga
gtcccagttactgccagatgtgggagctattgtaacctgtaagggtctctagcatcaatfcacgctttgccaaagtacacatctgtatgtgggtccatgcctcttaagaactcttttcgaggaaactccgcaagggaagatgtccgagcaactgaaaaagacaaggttgaaattataagagtttc
cgcccaggtgacattgtcttgccaaagtgtatctcttaggtgatgcacagtcgaactacgtgtaaccacgcccagagacgagctgggagtggtgtagccacagtgagtcaggtatccagatggttcccatcagctggtgtgagatgcagtgccctaagaccacactaaagaattcc
ggaaagtagcccgagtacaacccgaattctgcagacctaa