

KLK11-kallikrein-related peptidase 11 Gene

Reference: PTXBC022068

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRSS20; TLSP; kallikrein-11; hK11; hippostasin; serine protease 20; trypsin-like protease; kallikrein related peptidase 11		

For research use only.

GENE DESCRIPTION	kallikrein-related peptidase 11
GENE ID	11012
ORIGIN SPECIES	Human
ACCESSIONS	BC022068
VECTOR	pDONR223
NCBI SYMBOL	KLK11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggattctgcagttatcctgctgctgcaacagggttagggggagagaccaggatcatcaagggttcgagtgcaagcctcactccagccctggcaggcagccctgttcgagaagacgcggctactctgtggggcgacgctcatcgccccagatggctcctgacagc
agccctgctcctcaagccccgctacatagttcactggggcagcacaacctccagaaggaggagggctgtgagcagaccggacagccactgagtcctccccacccccggctcaacaacagcctcccaacaaagaccacccgcaatgacatcatgctgtgaagatggcatgcc
agtcctcatcactgggctgtgcgacccctcaccctctcctcagctgtgtgctactgctggcaccagctgctcatttcgggctggggcagcagctccagccccagttacgctgctcctcacaccttgcgatgcgccaacatcaccatcattgagcaccagaagtgtgagaacgctaccccg
gcaacatcacagacacatgggtgtgtgccagctgcagggaaggggcaaggactcctgccagggtgactccggggccctctggtctgtaaccagctcttcaaggcattatctctggggccagatccgtgtgcgatcacccgaaagcctgtgtctacacgaaagtctgcaaatatg
tggactgcatccaggagacgatgaagaacaattag